

		GLOVE MATERIAL	
	NITRILE	LATEX	
Acetaldehyde: acetic aldehyde	Not recommended	Good	
Acetic acid 50%	Fair	Excellent	
Acetic acid, glacial	Fair	Good	
Acetic anhydride	Fair	Good	
Acetone	Not recommended	Fair	
Acetonitrile	Fair	Fair	
Acrylic acid	Good	Good	
Alcoholic beverages	Excellent	Excellent	
Ammonium acetate	Excellent	Excellent	
Ammonium carbonate	Excellent	Excellent	
Ammonium chloride	Excellent	Excellent	
Ammonium concentrate	Good	Excellent	
Ammonium fluoride 30-70%	Excellent	Excellent	
Ammonium hydroxide <30%	Excellent	Excellent	
Ammonium hydroxide 30-70%	Excellent	Excellent	
Ammonium nitrate	Excellent	Excellent	
Amylic Alcohol	Good	Fair	
Aniline	Fair	Not recommended	
Animal Fats	Excellent	Fair	
Aqua Regia	Fair	Fair	
Asphalt	Excellent	Not recommended	
AZT	Not recommended	Good	
<b>B</b> eet	Excellent	Excellent	
Benzaldehyde: benzoic aldehyde	Fair	Not recommended	
Benzene	Fair	Not recommended	
Benzyl alcohol	Fair	Fair	
Bleach	Excellent	Good	
Borax	Excellent	Excellent	
Boric acid	Excellent	Good	
Brake fluid: lookheed	Excellent	Fair	
Bromides	Excellent	Fair	
Bromopropionic acid	Fair	Good	
Butoxyethanol	Excellent	Good	
Butter	Excellent	Not recommended	
Butyle acetate	Good	Not recommended	
Butyle cellusolve	Good	Good	
Calcium chloride	Excellent	Excellent	
Calcium disulfide	Good	Fair	

CHEMICAL	GLOVE MATERIAL	
	NITRILE	LATEX
Calcium hydroxide	Excellent	Excellent
Calcium hypochloride	Excellent	Excellent
Calcium nitrate	Excellent	Excellent
Calcium oxide	Excellent	Excellent
Calcium phosphate	Excellent	Excellent
Carbon tetrachloride	Good	Not recommended
Castor oil	Excellent	Not recommended
Chlorine	Excellent	Fair
Chloroacetone	Not recommended	Excellent
Chlorobenzene	Fair	Fair
Chlorodibromomethane	Fair	Fair
Chloroform	Fair	Not recommended
Chloronaphthalenes	Fair	Fair
Chromic acid	Fair	Fair
Cisplatin	Good	Good
Citric acid	Excellent	Excellent
Creosote	Excellent	Fair
Cresol	Excellent	Good
Cutting oil	Excellent	Not recommended
Cyclohexane	Excellent	Not recommended
Cyclohexanol	Excellent	Excellent
Cyclohexanone	Not recommended	Good
Cyclohexylamine	Fair	Fair
<b>D</b> i isobutyl ketone	Good	Fair
Di isobutylamine	Excellent	Fair
Di methyl ether	Good	Fair
Di methyl sulfoxide: DMSO	Good	Excellent
Di methylacetamide	Fair	Good
Di methylformamide: DMF	Fair	Fair
Diacetone alcohol	Good	Excellent
Diallylamine	Fair	Fair
Dibutyl phtalate	Excellent	Excellent
Dibutylether	Good	Not recommended
Dichloroacetyl chloride	Fair	Fair
Dichloroethane	Fair	Not recommended
Diesel oils	Excellent	Not recommended
Diethanolamine	Excellent	Excellent
Diethylamine	Good	Fair



	GLOVE MATERIAL	
	NITRILE	LATEX
Diethylene glycol	Excellent	Excellent
Diethylenetriamine	Fair	Fair
Di-n-amylamine	Excellent	Fair
Di-n-butyl phthalate	Excellent	Fair
Di-n-butylamine	Excellent	Fair
Di-n-octyl phthalate	Excellent	Fair
Dioctyl phtalate	Excellent	Fair
Dyes: hair	Excellent	Excellent
1, 3-Dioxane	Fair	Fair
1, 4-Dioxane	Fair	Fair
Epichlorohydrin	Fair	Fair
Ethanol: Ethyl Alcohol	Excellent	Good
2-Ethoxyethanol	Excellent	Fair
2-Ethoxyethylacetate	Fair	Not recommended
Ethyl acetate	Fair	Fair
Ethyl ether	Good	Fair
Ethylaniline	Excellent	Fair
Ethylene dichloride	Fair	Fair
Ethylene glycol	Excellent	Excellent
Ethylene glycol	Excellent	Excellent
Ethylene glycol dimethyl ether	Fair	Fair
Fertiliser	Excellent	Excellent
Fish & Shellfish	Excellent	Fair
Fixing agents	Excellent	Excellent
Flourides	Excellent	Fair
Formaldehyde 30% - 70%	Excellent	Good
Formic acid	Good	Excellent
Freon 113 or TF	Excellent	Fair
Freon TMC	Fair	Fair
Fuels	Excellent	Not recommended
Furaldehyde	Not recommended	Good
Furfura	Fair	Fair
Gas oils	Excellent	no
Gasoline, 40-50% aromatics	Excellent	Fair
Gasoline, unleaded	Good	Fair
Glutaraldehyde <5%	Good	Good
·	Excellent	Excellent
Slycerine Slycerophtalic paint		
Slycerophtalic paint	Excellent	Not recommended
Glycols	Excellent	Excellent

CHEMICAL	GLOVE MATERIAL	
	NITRILE	LATEX
Heptanes	Excellent	Fair
Hexamethyldisiloxane	Good	Fair
Hexane	Excellent	Not recommended
Household Detergents	Good	Excellent
Hydraulic fluid: Esters	Excellent	Excellent
Hydraulic fluid: Petrol	Excellent	Not recommended
Hydrazine	Excellent	Fair
Hydrochloric acid <30%	Good	Excellent
Hydrochloric acid 30-70%	Good	Good
Hydrofluoric acid <50%	Excellent	Excellent
Hydrogen peroxide	Excellent	Fair
Isobutanol: isobutylic alcohol	Excellent	Good
Isobutylcetone	Not recommended	Excellent
Isooctane	Excellent	Fair
Isopropyl alcohol	Excellent	Excellent
Isopropylamine	Good	Good
Jet fuel <30% aromatics 73-248C	Good	Good
Kerosene	Excellent	Not recommended
Lactic acid 85%	Good	Good
Lard oil	Excellent	Not recommended
Lauric acid	Good	Good
Linseed oil	Excellent	Not recommended
Lubricating oil	Excellent	Not recommended
Magnesium oxide	Excellent	Excellent
Malathion 30-70%	Excellent	Excellent
Maleic acid	Good	Good
Methanol: methyl alcohol	Excellent	Fair
2-Methoxyethanol	Excellent	Fair
Methyl acetate	Fair	Fair
Methyl ethyl ketone	Not recommended	Good
Methyl isobutyl ketone	Not recommended	Good
Methyl methacrylate	Fair	Fair
Methylamine	Excellent	Good
Methylaniline	Excellent	Fair
Methylene chloride	Fair	Not recommended
Milk & dairy products	Excellent	Fair
Mineral fats	Excellent	Not recommended
Monochlorobenzene	Not recommended	Fair
Monoethanolamine	Excellent	Excellent
Naphta: white spirit <3% aromatics	Excellent	Fair



	GLOVE MATERIAL	
	NITRILE	LATEX
Naphta: white spirit 15-20% aromatics	Excellent	Fair
Naphtalene	Good	Not recommended
n-Amyl acetate	Fair	Fair
n-Butanol: butylic alcohol	Excellent	Good
n-Butyl acetate	Fair	Fair
n-Butyl alcohol	Excellent	Excellent
Nitric acid <30%	Excellent	Excellent
Nitric acid 30-70%	Fair	Good
Nitrobenzene	Fair	Good
Nitroethane	Fair	Excellent
Nitrohydrochloric acid	Fair	Not recommended
n-Methyl-2-Pyrrolidone	Fair	Excellent
n-Nitrosodiethylamine	Fair	No Information
1-Nitropropane	Fair	Good
2-Nitropropane	Fair	Fair
Non alcoholic beverages	Excellent	Excellent
n-Propyl alcohol	Excellent	Excellent
Octane	Excellent	Excellent
Octanol: octyl alcohol	Excellent	Excellent
Oils for turbines	Excellent	Not recommended
Oleic acid	Excellent	Good
Olive oil	Excellent	Not recommended
Oxalic acid	Excellent	Excellent
Palmitic acid	Excellent	Fair
Paraffin oil	Excellent	Not recommended
Peanut oil	Excellent	Not recommended
Pentachlorophenol	Good	Fair
Pentane	Excellent	Fair
Perchlorethylene	Good	Fair
Perchloric acid 30-70%	Excellent	Fair
Perfumes & essences	Excellent	Excellent
Peroxyacetic acid	Fair	Fair
Petrol	Excellent	Not recommended
Petrolium ethers: 80-110C	Good	Fair
Petrolium products	Good	Not recommended
Phenol: phenic alcohol	Fair	Good
Phenol: phenic alcohol >70%	Excellent	Fair
Phosphoric acid 75%	Excellent	Good
Picric acid	Excellent	Good
Polychlorinated biphenyls: bcp	Good	Fair

	GLOVE MATERIAL	
	NITRILE	LATEX
Polyester resins	Good	Not recommended
Potassium bicarbonate	Excellent	Excellent
Potassium bichromate	Excellent	Fair
Potassium carbonate	Excellent	Excellent
Potassium carbonate conc.	Excellent	Good
Potassium chloride	Excellent	Excellent
Potassium cyanide	Excellent	Excellent
Potassium hydroxide	Excellent	Good
Potassium iodide	Excellent	Excellent
Potassium nitrate	Excellent	Excellent
Potassium permanganate	Excellent	Excellent
Potassium phosphate	Excellent	Excellent
Potassium sulphate	Excellent	Excellent
Poultry	Excellent	Fair
Propyl acetate	Fair	Fair
Pyridine	Fair	Fair
Setting agents	Excellent	Excellent
Shampoos	Excellent	Excellent
Silicate	Excellent	Excellent
Silicone etch	Fair	Fair
Silver nitrate	Good	Excellent
Sodium bicarbonate	Excellent	Excellent
Sodium carbonate	Excellent	Excellent
Sodium chloride	Excellent	Excellent
Sodium fluoride	Excellent	Excellent
Sodium hydroxide 30-70%	Excellent	Excellent
Sodium hypochloride	Excellent	Excellent
Sodium hypochlorite	Excellent	Excellent
Sodium nitrate	Excellent	Excellent
Sodium phosphates	Excellent	Excellent
Sodium sulphate	Excellent	Excellent
Sodium thiosulfate	Excellent	Excellent
Soldium bisulphate	Excellent	Excellent
Soya bean oil	Excellent	Not recommended
Styrene	Fair	Not recommended
Sulphites: bi and hypo	Excellent	Excellent
Sulphuric acid <30%	No Information	Excellent
Sulphuric acid >70%	Fair	Fair
Sulphuric acid 30-70%	Fair	Excellent
Tannic acid	Good	Good



	GLOVE MATERIAL	
	NITRILE	LATEX
THF: Tetrahydrofurane	Fair	Not recommended
Toluene	Good	Not recommended
Toluene-2,4-Diisocyanate: TDI	Fair	Fair
Tributylphosphate	Not recommended	Not recommended
Trichlorethylene	Fair	Not recommended
Tricresyl phosphate	Good	Good
Triethanolamine 85%	Excellent	Excellent
Trinitrobenzine	Good	Not recommended
Trinitrotoluene	Good	Not recommended
Triphenylphosphate	Fair	Not recommended
Turnipseed oil	Excellent	Not recommended
Turpentine	Excellent	Not recommended
Turpentine spirit	Excellent	Not recommended
1,2,4,5- Tetrachlorobenzene	Excellent	No Information
1,1,1,2- Tetrachloroethane	Fair	Fair
1,2,4-Trichlorobenzene	Fair	Fair
1,1,1-Trichloroethane	Fair	Fair
1,1,2-Trichloroethane	Fair	Fair
Vinegar & condiments	Excellent	Excellent
Vinyl acetate	Fair	Not recommended
<b>W</b> ashing powders	Excellent	Excellent
Water paints	Excellent	Excellent
Weedkillers	Excellent	Good
<b>X</b> ylene	Good	Not recommended
Xylophene	Good	Not recommended
<b>Z</b> inc sulphate	Excellent	Excellent

Understanding the use of Thin Wall gloves: TechNiGlove products are not intended for immersion or prolonged contact with harmful substances, but are manufactured for use against incidental contact or "splash protection". This Chemical Resistance Chart is intended as a general guide, it is the responsibility of the user to choose the appropriate type of glove and thickness for their individual application. This Guide is based on public domain data and TechNiGlove products have not been individually tested or verified for use against the chemicals listed.

Please contact our Customer Service Team to choose the right Techniglove for your specific application.